

PROTOCOL FOR SEQUENCING PCR OF SARS-COV-2

This is a protocol for multiplex PCR reactions prior to the library preparation of the Oxford Nanopore sequencing.

Materials

- o primers pool A and B
- o Micropipettes and tips
- o Microcentrifuge tubes size 1,5mL
- o PCR microtubes with attached caps size 0.2 mL
- o Disposable lab coats
- o Disposable powder-free lab gloves
- o 0,2mL microtube racks
- o Thermocycler

Protocol version
updated on 01 June 2021

Mastermix

All primers pools should be diluted to a total of 10uM from the 100uM stock.

ARTIC NETWORK PROTOCOL – Quick, Josh, 2020

Components	Volume (uL)		Components	Reduced Volume (uL)
Nuclease Free water	12,75		Nuclease Free water	5,33
Q5 reaction buffer (5x)	5		Q5 reaction buffer (5x)	2,5
dNTPs (10mM)	0,5		dNTPs (10mM)	0,25
Primers for the pool 10uM	4		Primers for the pool 10uM	1,8
Q5 DNA polymerase	0,25		Q5 DNA polymerase	0,125

Prepare **2** mastermixes with the above components for the appropriate number of samples. In one tube add the **primers** of the pool **A** and in the other tube, add primers of the pool **B**.

After the mastermix is finished, pipette out **22.5uL (or 10 uL if using the reduced reaction)** of the **mastermix** into PCR microtubes, and then add in each tube **2.5ul** of the appropriate cDNA.

Thermocycling parameters (total time: ~2h37m)

Step	Temperature	Time	Cycle number
Initial denature	98°C	30 seconds	1
Denature	98°C	15 seconds	30
Annealing and extension	65°C	5 minutes	
Hold	4°C	∞	1